
Phosphate Buffer Solution for Giemsa Staining (PFG)

In vitro diagnostic medical aid.
Only for professional use.

Purpose: The product is designed for the staining of dry coatings, imprints, blood smears.

Package contents: 500 ml solution in a bottle

Composition: disodium hydrogen phosphate
sodium dihydrogenphosphate
demineralized water

Storage: Store in a cool place at +2 to +8 °C.

Usability: Use within 30 days after opening the primary packaging.

Procedure: Use the solution according to internal procedures.

Information for users: Do not use after the expiration date.

Dispose of used and unused products as hazardous waste in accordance with applicable regulations.

Instructions in case of a damaged packaging: In case of damage to the packaging, do not use and dispose of the product.

Manufacturer

EnviroLab s.r.o.
Galvaniho 7
821 04 Bratislava
Slovak Republic

Place of business:

Obchodná ulica 7
946 51 Nesvady
Slovak Republic